

Veterinary Lasers from Excel Lasers



AngioDynamics Delta laser

Excel Lasers is pleased to be able to offer a limited quantity of reconditioned AngioDynamics Delta surgical lasers, exclusively to veterinary surgeons. These medical-grade, 810 nm diode laser systems are available with power outputs of 15 W or 30 W. They have been fully serviced and calibrated, and are offered with a full 12 months warranty.

These portable desktop units, designed and manufactured in the United Kingdom, are ideal for incision, ablation and coagulation of soft tissue. Energy is delivered through optical fibres which are functional in fluid environments and by application through an endoscope.

- Fully-reconditioned & serviced
- 12 months Excel Lasers warranty
- Foot switch & accessory pack
- Designed and built in the UK.
- Installation and training available
- Full range of fibres and accessories

Specifications	DELTA 15	DELTA 30
Maximum Power	15 W	30 W
Wavelength	810 nm	
Operational Modes	CW, Single or Repeat Pulse	
Pulse Duration	100 ms – 9900 ms	
Pulse Interval	100 ms – 1000 ms	
Fibres (SMA905)	400 µm, 600 µm or 1000 µm (>0.37 Numerical Aperture)	
Aiming beam	Red Class 3R diode laser, ≤ 5 mW, 635-660 nm	
Voltage/Current Rating	100 – 240 V AC, 50 – 60 Hz, 500 VA max.	
Dimensions	370(W) x 325(D) x 150(H) mm	
Weight	12 kg	

Distributed and supported in the UK by:

Excel Lasers Limited
 The Old Yard
 Great Bradley Hall
 Newmarket, CB8 9LT
 United Kingdom

EXCEL LASERS 

Tel: +44 (0)1440 783201
 Email: sales@excellasers.com
 WWW: www.excellasers.com